MD 555 cleaner

Special detergent for suction systems



Foam-free special cleaner for dental suction units including drainage lines for regions with problematic tap water and drainage lines with insufficient gradient.



Contact time Recommendation

5 %
30-120
min
at least twice weekly

dissolves inorganic deposits of prophy powder and pearl products

maintains the value of your suction system

tested material compatibility

recommended by leading equipment manufacturers

intensive cleaning effect

works perfectly in combination with Orotol plus

Function



"Descaling agent/intensive cleaning agent": excellent cleaning performance to tackle deposits, including limescale resulting from very hard water, stubborn salts, etc.



Key features

5%

application concentration

intensive cleaning
effect to remove
deposits and
accumulations of silt

guarantees
technically and
hygienically
perfect operation of
the suction system

